

ISO 9001 : 2000 Company

**NULIFE**

COMPRESSOR ▪ NEBULIZER FOR  
**Aerosoltherapy**

CE  
0434

**SLIMNEB<sup>TM</sup>**

SPACE SAVING WITH ACCESSORIES COMPARTMENT AND RECESSIVE HANDLE

### PISTON TYPE

- Low Noise Level
- Safe for Long Use
- Recessive Handle
- Accessories Compartment



**NULIFE** and **YOU**  
TOGETHER LET'S SAVE LIVES



# SLIMNEB<sup>®</sup>

SPACE SAVING WITH ACCESSORIES COMPARTMENT AND RECESSIVE HANDLE

## Breath easy with range nebulizers from **NULIFE**

With **SLIMNEB**<sup>®</sup> nebulizers, an effective treatment of respiratory diseases can be achieved. We have a wide range of suitable devices & accessories for every individual need.

Every **SLIMNEB**<sup>®</sup> nebulizer is made for long lasting use while maintaining the highest standard of quality. A high aerosol output and well-controlled nebulization of precious medication make our products popular amongst many private and professional users.

### **SLIMNEB**<sup>®</sup>

- Very powerful & reliable piston compressor nebulizer especially for professional & intensive use.
- Powerful compressor and matched nebulizer kit, ensures medication is atomized into fine particles faster, reaching respiratory tracts more effectively and resulting in a better medication compliance.
- Equipped with safety fuse for safety.
- Wide Range of accessories available.
- Ample storage space for accessories.
- Anti-shock case with recessive handle for easy carrying.

- ❖ Type of compressor ..... Piston
- ❖ Pressure ..... Max 2.5 bar
- ❖ Air Flow ..... Min 8 LPM
- ❖ Noise at 1 meter ..... 55 dBA
- ❖ Operating times ..... 60 min ON 30 min OFF
- ❖ Particle sizes ..... 0.5 - 5 microns
- ❖ Size ..... 140 mm x 104 mm x 290 mm
- ❖ Power Supply ..... (215 - 240vAc; 50 - 60Hz)  
(110 - 130vAc, 50-60Hz)
- ❖ Power Consumption .... 90 watt maximum
- ❖ Weight ..... approx 1.6 kg

\* Accessories supplied with the equipment may vary.

Specification : Subject to change without prior notice.



Wide range of accessories available\*



### **MRK HEALTHCARE**

An ISO 9001 : 2000 Certified Company

B4/ B5, Byculla Service Industries, D.K. Marg,  
Byculla, Mumbai 400 027. INDIA.

Tel. : 91-22-237 48371/2/3 Fax : 91-22-2374 8370

E-mail : [info@mrkhealthcare.com](mailto:info@mrkhealthcare.com)

Website : [www.mrkhealthcare.com](http://www.mrkhealthcare.com)

